

Aim :-

To determine surface tension of given liquids by drop count method.

Reference :Requirement :->

- (a) Apparatus :- Pycnometer, Stalagmometer, Pipette, Thermometer, Analytical balance, Weight box, beaker, etc.
- (b) Chemicals :-> Distilled water.

Theory :

Surface tension: This property of liquids arises from the intermolecular forces of attraction. A molecule in the interior of a liquid is attracted equally in all directions by the molecules around it.

Its SI unit is newton per metre (Nm^{-1}).

Teacher's Signature _____

Surface Tension Determination methods :

The methods commonly employed for the determination of surface tension are :-

1. Capillary-rise method
 - (a) Single capillary-rise method
 - (b) Double capillary-rise method
2. Drop formation method
 - (a) Drop ~~no~~ number method
 - (b) Drop weight method
3. Ring-detachment method
4. Maximum bubble pressure method
5. Wilhelmy plate method.
6. Pendant drop method.

Procedure:

1. Thoroughly clean pycnometer and stalagmometer with chromic acid and wash two times with fresh distilled water.
2. Stalagmometer must be fixed in a vertical position using stand.

Teacher's Signature _____

Observations :

1. Room temperature = °C
2. Weight of Empty pycnometer = W_1 .
3. Weight of Pycnometer + Distilled water = W_2
4. Weight of Pycnometer + Liquid = W_3 .

Liquids	Number of drops				Specific gravity	Surface tension.
	I	II	III	Mean		

Calculations :

1. Weight of liquid = $W_3 - W_1$.
2. Weight of distilled water = $W_2 - W_1$.
3. Specific gravity of liquid = $\frac{(W_3 - W_1)}{(W_2 - W_1)}$

$$\gamma_2 = \frac{P_2 n_2}{P_1 n_1}$$

Calculate surface tension of other liquids by substituting data in place of distilled water using same equation.

3. Fill water in stalagmometer up to mark A and count total number of drops formed from mark A to B.
4. Repeat step 3 at least 3 times for accuracy.
5. Wash stalagmometer using same liquid of which surface tension is to be determined.
6. Repeat step 3, 4, 5 for other chemicals.
7. Determine density of liquids by using pycnometer.

Result:

Surface tension of the given liquid has been determined in the laboratory successfully.